Speed Duo Diro/Leish K™

Speed Duo Diro/Leish K™

Veterinary diagnostic kit

For *in vitro* use only

■ CLINICAL BENEFIT

Heartworm infection is a helminthiasis that affects dogs and cats, and is caused by the presence and development of a filaria called *Dirofilaria immitis* in the pulmonary arteries and right ventricle. It is transmitted by mosquito bites. Heartworm infection causes a cardiac cough, heart murmur, exercise intolerance, and dyspnea. Testing for circulating antigens reveals the presence of adult worms; occult *Dirofilaria* can be found in the absence of circulating microfilaria (30% of cases).

Canine leishmaniosis is caused by the proliferation of a parasite, *Leishmania infantum*, in the macrophages of its host. It is transmitted by a small dipter: the sandfly (*Phlebotomus spp.*). Leishmaniosis is responsible for a range of clinical signs due to an immune reaction with variable consequences on the organs: depression, lymphadenopathy, anaemia, proteinuria, etc.

Heartworm infection and Leishmaniosis are found almost worldwide, associated with the distribution of their vector. In Europe, both diseases are particularly frequent around the Mediterranean basin. After infection, Heartworm infection and Leishmaniosis present with relatively long asymptomatic or subclinical phases, making their diagnosis difficult. In endemic areas, given the potential co-existence and severity of these diseases, it is advisable to perform regular, systematic screening to adapt preventive measures.

■ PRINCIPLE

Speed Duo Diro/Leish K is a rapid qualitative test, using membrane-based immunochromatography, enabling the detection of circulating antigens of *Dirofilaria immitis* and anti-*Leishmania infantum* antibodies in dogs.

For each test, one drop of sample (serum, plasma, or anticoagulated whole blood) is added to the sample well. In the sample well, the stained particles of the conjugate bind to any *Dirofilaria immitis* antigens (Diro part) and to any anti-*Leishmania infantum* antibodies (Leish part) present in the sample. The resulting conjugate/anti-*Leishmania infantum* antibody and conjugate/*Dirofilaria immitis* antigen complexes migrate along the membrane via capillarity. They are captured by specific anti-*Dirofilaria immitis* antibodies and specific *Leishmania infantum* antigens that are bound to the membrane, forming a pink test band by accumulation of coloured particles (central part of the window). The mixture continues to migrate to the end of the membrane, where the remaining coloured particles form a pink control band, confirming the validity of the test.

■ OPERATING PROCEDURE

► FOR EACH TEST YOU WILL NEED:

1 test kit, 1 single-use pipette, and the bottle of reagent.

Use all reagents at room temperature.

Never mix reagents from different batches.

1/ ADD THE SAMPLE:

-Serum, plasma, or anticoagulated whole blood (EDTA or heparin):

using a single-use pipette held vertically, **add 1 drop** of sample to each sample well.

2/ ADD THE REAGENT:

- Add the reagent immediately after the sample. Hold the bottle of reagent vertically and **add 5 drops** of reagent to each sample well.

- If no migration is observed within 2 minutes, **add an extra 2 drops** of reagent to the sample well.

3/ READING AND INTERPRETING THE RESULTS:

Wait 15 minutes for migration and read:

- **A NEGATIVE TEST:** any test showing only **1 pink band** (control band) in the reading window is a negative result.

- **A POSITIVE TEST:** any test showing **2 distinct pink bands** (test band + control band) is a positive result.

The apparition of a test band within only 10 minutes of migration is indicative of positive test. Any colour change in the test band, even faint, should be considered as a positive result.

- **The absence of control band makes the test invalid.**

■ RECOMMENDATIONS

• SHELF LIFE:

- 24 months at room temperature from the manufacturing date. The expiry date is indicated on the kit and each sachet of test module.
- Store at room temperature, between +2°C and +30°C. Avoid exposing the test to high temperatures or temperatures below freezing.

• SAMPLES :

- The samples of serum, plasma, or anticoagulated whole blood should be at room temperature.
- Storage of samples for use within 24 hours: store serum, plasma, or anticoagulated whole blood at room temperature.
- Storage of samples for use within 72 hours: store serum or plasma at between +2°C and +8°C (refrigerator). The storage of anticoagulated whole blood for several hours between +2°C and

+8°C can lead to partial haemolysis (not visible to the naked eye), which may interfere with the interpretation of the results.

- Testing more than 72 hours after sampling: store serum and plasma at -20°C. Never freeze a whole blood sample.

• **PRECAUTIONS OF USE :**

- To ensure correct migration of the sample, place the test device on a flat, horizontal surface.

- To standardise the amount of sample recommended for each test, only use the pipettes supplied in the kit.

- To ensure that the sample and reagent are correctly applied, hold the reagent dropper bottle and pipette vertically.

• **OTHER RECOMMENDATIONS :**

- When performing the test, the reading window may show a slight pink background coloration; this does not affect the quality of the results.

- **Do not mix reagents from different batches.**

- Allow the reagents to come to room temperature before use.

- Use a new single-use pipette for each test.

The above recommendations are only guidelines; no test is 100% accurate at all times and under all conditions. The purpose of this test is to detect *Dirofilaria immitis* antigens and *Leishmania infantum* antibodies to diagnose Heartworm infection and Leishmaniosis respectively. All test results should be interpreted in the light of the patient's history, clinical examination, and the results of any further diagnostic tests. The definitive diagnosis remains the veterinarian's prerogative and responsibility.

Bio Veto Test and its distributors cannot be held responsible for the consequences of misuse or misinterpretation of the results given by the test.

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