**Speed FIV™**

# Speed FIV™

## Veterinary diagnostic kit

### For *in vitro* use only

#### ■ CLINICAL APPLICATION

Feline Immunodeficiency virus (FIV) is a retrovirus that is present worldwide. It is primarily transmitted through close contact between cats via saliva or bites.

The initial viraemia stabilises within a few weeks of infection with FIV, stimulating the production of specific antibodies directed against the viral proteins. The circulating viral load decreases; this is followed by an asymptomatic phase, which can last several months to several years. The animal is then considered as a healthy carrier, with a potential for viral excretion.

During the clinical phase of FIV infection, symptoms are polymorphic and non-specific: depression, anorexia, pyrexia, immunodeficiency syndrome, anaemia, etc. Further diagnostic tests are required to obtain a definitive diagnosis.

Regular FIV testing can be used to determine disease status and to indicate whether specific precautions need to be taken (separation from healthy cats, neutering, and medical supervision...).

#### ■ PRINCIPLE

Speed FIV is a rapid qualitative immunoassay, using immunochromatography, for the detection of anti-FIV antibodies in cats.

For each test, one drop of the sample (serum, plasma or

anticoagulated whole blood) is added to the sample well. In the sample well, the stained particles of the conjugate bind to any anti-FIV antibodies present in the sample. The resulting conjugate/antibody complexes migrate along the membrane via capillarity. They are captured by specific FIV antigens that are bound to the membrane, forming a pink test band by accumulation of coloured particles. The mixture continues to migrate along the strip to the end of the membrane where the remaining coloured particles form a pink control band, which confirms the validity of the test.

## ■ PROCEDURE

### ► FOR EACH TEST YOU WILL NEED:

1 test device, 1 single-use pipette and the reagent dropper bottle. The reagents should be used at room temperature.

**Do not mix reagents from different batches.**

### 1/ ADD THE SAMPLE:

Serum, plasma, anticoagulated whole blood (EDTA or heparin): with a single-use pipette held vertically, **add 1 drop** of sample to the sample well.

### 2/ ADD THE REAGENT:

- Add the reagent immediately after the sample. Hold the dropper bottle in a vertical position and **add 5 drops** of reagent to the sample well.

- If no migration is observed within 2 minutes, **add an extra 2 drops** of reagent to the sample well.

### 3/ READING AND INTERPRETATION OF RESULTS:

**Wait 15 minutes** for migration and read:

- **NEGATIVE TEST:** Any test showing only **1 pink band** (control band) is a negative result.

- **POSITIVE TEST:** Any test showing **2 distinct pink bands** (test band + control band) indicates a positive result.

The test is positive if a pink test band appears within 10 minutes of migration. Any colour change in the test band, even faint, should be considered as a positive result.

- **If no control band appears, the test is invalid.**

### ■ RECOMMENDATIONS

#### • SHELF LIFE:

- 24 months at room temperature (from manufacturing date). The expiry date is printed on each box and on each test pouch.
- Store at room temperature, between +2°C to +30°C. Keep the test away from sources of excessive heat or cold.

#### • SAMPLES:

- Serum, plasma, or anticoagulated whole blood must be used at room temperature.
- For testing within 24 hours of sampling: store serum, plasma, or anticoagulated whole blood at room temperature.
- For testing between 24 to 72 hours of sampling: store serum or

plasma between +2°C to +8°C (refrigerator). Storing whole blood between +2°C to +8°C for several hours may induce haemolysis (not visible to the naked eye), which may interfere with the interpretation of the results.

- Testing more than 72 hours after sampling: store serum or plasma at -20°C (freezer). Never freeze a whole blood sample.

• **PRECAUTIONS:**

- To ensure correct migration of the sample, place the test device on a horizontal surface.

- To standardize the amount of sample recommended for each test, only use the single-use pipettes supplied in the kit.

- To ensure that the sample and reagent are correctly applied, always hold the pipette or dropper bottle vertically.

• **OTHER RECOMMENDATIONS:**

- The reading window may show a light pink background colouration; this does not affect the quality of the results.

- **Do not mix reagents from different batches.**

- Allow the reagents to come to room temperature before use.

- Use a new single-use pipette for each test.

The above recommendations are only guidelines; no test is 100% accurate at all times and under all conditions. The purpose of this test is to detect anti-FIV antibodies to diagnose FIV. All test results should be interpreted in the light of the patient's history, clinical examination, and the results of any further diagnostic tests. The definitive diagnosis remains the veterinarian's prerogative and responsibility.

Bio Veto Test and its distributors cannot be held responsible for the consequences of misuse or misinterpretation of the results of the test.

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