Speed V-Diar 5TM

English	1
Français	5
Español	9
Português	13
Italiano	17
Ελληνικά	21
Deutsch	25
Netherlands	29
Procedure	33

FOR IN VITRO USE ONLY ENGLISH

■ CLINICAL INTEREST

Neonatal calf diarrhoea is a common problem that primarily arises during the first three weeks of life. The resulting reduction in growth rates and the rapid contamination of the farm during the calving season are major economic concerns for the farmer. However, although the cause cannot be determined by symptoms alone, determining an etiological diagnosis right from the start of an outbreak is essential to direct therapeutic and prophylactic measures.

Speed V-Diar 5 is a rapid test that enables the veterinarian to simultaneously test for the five main agents responsible for neonatal calf diarrhoea: Coronavirus, Rotavirus, *Cryptosporidium parvum*, E. coli F5 (K99) and E. coli CS31A.

■ PRINCIPLE

Speed V-Diar 5 is a rapid immunochromatographic strip test, for the detection of specific antigens of the five main agents responsible for neonatal diarrhoea in the calf: Coronavirus, Rotavirus, *Cryptosporidium parvum*, E. coli F5 (K99) and E. coli CS31A. There is one test strip for each of these pathogenic agents. The five test strips are grouped together in one device enabling combined testing.

For each test strip, the coloured particles of the conjugate, integrated into the membrane, bind to the antigens of any of the five pathogens present in the sample. The resulting conjugate/ antigen complexes migrate via capillarity. They are then captured by the specific antibodies bound to the membrane, forming a red-coloured test band through the accumulation of coloured particles. The mixture continues to migrate to the end of the strip, where the remaining coloured particles form a red control band, which confirms the validity of the test.

■ OPERATING PROCEDURE

► FOR EACH TEST YOU WILL NEED:

One Speed V-Diar testing device containing 5 strips, one sample collection tube with reagent + sampling spoon, and one single use pipette.

Always use the reagents at room temperature.

Never mix reagents from different batches.

1/ PREPARING THE SAMPLE:

- Open the sampling tube containing the reagent.
- Place one volume of faeces in the reagent tube:
 - *For liquid faeces*: using the pipette, sample the faeces to be analysed and place **5 drops of the faeces** in the reagent tube.
 - *For non-liquid faeces*: using the spoon provided in the lid of the reagent tube, sample **one level spoonful of faeces** to be analysed and place it in the tube.
- Close the reagent tube again and mix the contents by inverting it a dozen times.

2/ ACTIVATING THE TESTING DEVICE:

- Unscrew the blue cap of the testing device.
- **Insert the tube of reagent in the testing device**, with the lid of the reagent tube uppermost.
- **Close the testing device** again using the blue lid, **screw shut until you hear a “click”**, and place it upright on a flat surface. Complete tightening of the lid of the testing device perforates the base of the reagent tube and allows its content to flow onto the device. On contact with the base of the test strips, the Reagent/Sample mixture will migrate vertically via capillarity on each test strip.
- If the liquid in the reagent tube does not start flowing after a few seconds, unscrew and re-screw the lid of the testing device slightly.
- **Leave it to migrate in the vertical position.**

3/ READING AND INTERPRETING THE RESULTS:

Wait **15 minutes for migration** and read.

The identification of the pathogenic agent is printed on the label at the top of each strip.

- **A NEGATIVE TEST** gives a single red band in the upper part of the strip (control band).

- **A POSITIVE TEST** gives two red bands (test band and control band).

The apparition of a test band after only 10 minutes of migration is evidence of a positive result for the valence in question. Any colouration of the test band, however slight, should be considered as a positive result.

- The absence of a control band in the upper part of the strip makes the test invalid.

■ RECOMMENDATIONS

• STABILITY / STORAGE:

- 18 months at room temperature from the date of manufacture. The expiry date is indicated on the box and on each testing device pouch.

- Store at room temperature, between +2°C and +30°C. Avoid exposing the test to high temperatures or below 0°C.

- If storing between +2°C and +8°C, allow the reagents to reach room temperature before use.

- The strips are sensitive to humidity and light. Only open the Speed V-Diar sachet when you are ready to test.

- Before use, the presence of visible green bands on the strips is normal.

- **SAMPLES:**

- It is preferable to perform the analysis on **fresh faecal matter**.
- If the test is delayed, store the faecal samples at between +2°C and +8°C and use within 24 hours. Beyond that, the faeces should be stored at -20°C.

- **PRECAUTIONS OF USE:**

- To standardise the quantity of sample needed for each test, only use the pipette or spoon provided with the kit.
- Do not take the test strips out of the testing device.
- It is advisable to wear disposable gloves when handling the samples and performing the test.

- **OTHER RECOMMENDATIONS:**

- **Never use reagent tubes and test devices from different kits.**
- Use a new test device, new tube of reagent, and a new sampling pipette for each test.
- Several pathogenic agents may be found in the same sample. Thus when a strongly positive result is seen on one of the strips, take the time to carefully examine the result on the other strips.

These recommendations are for guidance only, no diagnostic method is ever precise 100% of the time. The aim of this test is to help the veterinarian to make an etiological diagnosis of neonatal calf diarrhoea through the detection of soluble antigens of Coronavirus, Rotavirus, *Cryptosporidium parvum*, E. coli F5 (K99) and E. coli CS31A. When interpreting the test, the veterinarian should always consider the medical history and clinical and epidemiological background of the animal. The final diagnosis is the prerogative and responsibility of the veterinarian performing the test.

Bio Véto Test and its distributors cannot be held responsible for the consequences of incorrect use or interpretation of the results given by this test.

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