Fluoleish™

English	1
Français	5
Español	9
Português	13
Italiano	17
Ελληνικά	21
Deutsch	25
Netherlands	29

FOR *IN VITRO* USE ONLY ENGLISH

■ BENEFIT

Canine leishmaniosis is caused by the proliferation of a flagellated protozoan parasite, *Leishmania infantum*, in the macrophages of its host. The disease is transmitted by biting insects: sand flies, which serve as both host and vector.

The clinical manifestations of canine leishmaniosis are intermittent and of variable severity. This is why laboratory tests are needed to confirm the diagnosis. Serology is considered to be the most reliable diagnostic test as it provides an indication of the animal's immune response following an infestation. Several different serological tests are available; indirect immunofluorescence is considered the gold standard: it is a highly sensitive test and provides a quantitative evaluation of the animal's antibody levels. The intensity of the immunofluorescence reaction therefore increases during the phase of clinical expression, and monitoring over the time provides an indication of the progression of the disease.

■ PRINCIPLE

Fluoleish™ slides are used to perform an indirect immunofluorescence test to detect anti-*Leishmania* antibodies in the serum or plasma of the dog.

After depositing the sample to be analysed on one well sensitised by *Leishmania infantum*, any anti-*Leishmania* antibodies present in the sample will be captured by the leishmanias bound to the slide, and subsequently revealed by a canine anti-immunoglobulin that has been labelled with fluorescein.

The detection of antibodies is a qualitative test, any positive sample can then be quantitatively titrated using a simple series of dilutions.

■ MATERIAL REQUIRED

► Reagents needed to perform a Fluoleish™ test :

- one bottle of canine anti-immunoglobulin conjugated with fluorescein,
- one positive control,
- one negative control,
- PBS (Phosphate Buffered Saline),
- Tween® 80,
- Evans blue,
- Distilled water.

► **Required apparatus:** coverslips, an assembly medium for immunofluorescence, a humid chamber at +37°C, and a fluorescence microscope in UV light.

■ NATURE OF THE SAMPLE

Sample of canine serum or plasma.

■ OPERATING PROCEDURE

The reagents should be brought to room temperature 20 minutes before the start of the analysis.

► **For each test you will need:** 1 positive control well and 1 negative control well per slide, and as many wells as there are dilutions to be tested.

DILUTIONS OF THE SAMPLE:

- 1/50 Dilution: dilute 10 µL of the serum or plasma to be tested in 490 µL of PBS.
- 1/100 dilution and above: successive dilutions of 2 in 2 in PBS.

PERFORMING THE TEST:

1/ Place the following in separate wells:

- 15 µL of positive control,
- 15 µL of negative control,
- 15 µL of each dilution of the sample.

Allow to incubate for 30 minutes at +37°C in the humid chamber.

2/ Wash the slides by immersing them in a solution of PBS with Tween® 80: 2 drops of Tween® 80 for 1 litre of PBS (1 drop = 50 µL). Leave the slides immersed in the wash solution for 5 minutes under agitation.

Perform the second wash under agitation using a fresh solution for a further 5 minutes.

Briefly rinse in a bath of distilled water. Drain. Dry.

3/ In each well, place 15 µL of canine fluorescent anti-immunoglobulin conjugate diluted in a solution of PBS with Evans blue (1 drop - around 50 µL - of a 1% solution of Evans blue in 5 mL of PBS). The dilution of the conjugate should be defined for each batch of slides, from standard samples, as a function of the microscope used and the desired intensity of fluorescence.

Allow to incubate for 30 minutes at +37°C in the humid chamber.

4/ Wash the slides by immersing them in a solution of PBS with Tween® 80.

Leave the slides immersed in the wash solution for 5 minutes under agitation.

Perform a second wash under agitation using fresh solution for a further 5 minutes.

Briefly rinse in a bath of distilled water. Drain. Dry.

5/ Place 2 to 3 drops of the immunofluorescence mount medium on the surface of the slide and cover with a coverslip.

Read under the fluorescence microscope at 40 x 10.

■ READING

Negative reaction: absence of membrane fluorescence of leishmanias, which are stained dark red and are not easy to see.

Positive reaction: the leishmanias give off marked fluorescence, predominantly at the membrane. The flagellum is also fluorescent.

Intensity of the fluorescence of the sample tested: it can be more or less intense than the positive control.

Specific cases: in the event of a doubtful positive or negative control, it is important to repeat the procedure to validate your test samples under the right conditions.

■ INTERPRETATION

The titre of the sample corresponds to the last dilution that still shows fluorescence. The threshold of positivity is determined as a function of the standards and material used.

■ RECOMMENDATIONS

•STABILITY / STORAGE:

24 months between +2°C and +8°C.

•SAMPLES:

- Use samples of plasma or serum. It is preferable to work on fresh samples to prevent any risk of contamination or degradation.
- Storage for less than 5 days: store plasma and serum in the refrigerator between +2°C and +8°C.
- Storage for several weeks: store plasma and serum in the freezer at -20°C.

•GUIDELINES AND PRECAUTIONS OF USE:

- Do not use reagents that are past their expiry date.
- **Allow the slides and reagents to come to room temperature 20 minutes before the start of the analysis.**
- It is advisable to always perform 2 dilutions (1/50 and 1/100) to correctly evaluate the degree of basal fluorescence.
- It is advisable to perform one positive and one negative control well for each slide.
- If all of the wells of a slide have not been used for a test, the slide cannot be stored for future use.
- The dilution of the conjugate should be redefined, from standard samples, for each new batch of Fluoleish™ slides, as a function of the microscope used and the desired intensity of fluorescence.

These recommendations are for guidance only; no diagnostic method is ever precise 100% of the time. As with any diagnostic test, the results of the test should be interpreted in light of the history and clinical presentation.

Bio Véto Test cannot be held responsible for the consequences of incorrect use or interpretation of the results given by this test.

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