Speed ReSpiVB™

English	1
Français	5
Español	9
Português	13
Italiano	17
Ελληνικά	21
Deutsch	25
Netherlands	29
Procedure	33

FOR IN VITRO USE ONLY

ENGLISH

■ CLINICAL INTEREST

The frequency of respiratory problems in young cattle is high. However, the determination of the etiological diagnosis can sometimes be tricky: symptoms do not allow identification of the pathogen involved, which is essential for the establishment of appropriate therapeutic and prophylactic measures.

The Bovine Respiratory Syncytial Virus (BRSV) is the main pathogenic respiratory virus found in cattle under the age of 2 years. It is a contagious virus that can spread quickly in a herd if it is not identified early. The economic and medical consequences of respiratory disease in young bovines are severe: respiratory distress, reduced growth rates and high mortality rates (close to 25% for the younger animals).

BRSV is an intracellular, respiratory virus, which requires a nasopharyngeal sample containing many cells, and an extraction phase, to allow detection by immunological methods.

■ PRINCIPLE

Speed ReSpiVB is a rapid immunochromatographic strip test, for the detection of specific antigens of BRSV present in the respiratory epithelial cells of young cattle.

For each test, a deep nasopharyngeal sample is taken with a swab. The swab is discharged into a tube where an extraction step occurs. Then, a test strip is dipped into the tube containing the sample. The coloured particles of the conjugate bind to any antigens of BRSV present in the sample. The resulting conjugate/antigen complexes migrate by capillary action along the strip. They are then captured by the specific antibodies bound to the membrane, forming a pink-coloured test band due to the accumulation of coloured particles. The mixture continues to migrate to the end of the strip, where the remaining coloured particles form a pink control band, which confirms the validity of the test.

■ OPERATING PROCEDURE

► FOR EACH TEST YOU WILL NEED:

1 test strip, 1 test-tube, 1 swab, the bottle of distilled water, the bottle of extraction solution and 1 holder.

Always use the reagents at room temperature.

Never mix reagents from different batches.

1/ SAMPLING:

- A **deep nasopharyngeal swab** (or cytobrush).
- It is possible to perform the test on an individual sample (one swab per calf) or pooled samples from different calves (3 samples maximum) by discharging each swab in the same collection tube.

2/ PERFORMING THE SAMPLE:

- **Fill a sample tube with 5 drops of distilled water** for individual analysis (6 drops for a pooled analysis). Close the tube with the cap and put it in the holder.
- If you want to use a liquid obtained after trans-tracheal aspiration or bronchoalveolar lavage, quarter-fill a tube directly with the liquid harvested (250µL), without distilled water.
- Maintain the calf's head in extension, with nostrils raised. If they are obstructed with mucus, gently clean the openings of the nostrils with a compress.
- Remove the sterile swab from its protective packaging, taking care to conserve its sterility. **Introduce the whole length of the swab into one nostril, directing it upward to the dorsal meatus.**
- **Turn the swab and vigorously rub the mucous membrane** for several seconds until blood-tinged.
- Carefully remove the swab from the dorsal meatus and check that it is well saturated with "pink" secretions.
- In the minutes following the sampling, **insert the swab into the tube** containing the distilled water. **Discharge the swab** by crushing the tip against the inner sides of the tube. The swab can then be discarded. **Keep the tube containing the solution closed in order to perform the test in a clinic, within 2 hours (maximum).**

3/ SAMPLE EXTRACTION AND PERFORMING THE TEST:

- Add **4 drops from the extraction solution bottle** to the tube and close it. Gently homogenize and let the tube rest in the holder at room temperature for **10 minutes**.
- **Take one strip** from the bag containing the strips, and close it hermetically.
- Insert a strip into the tube, following the direction of the arrow.
- **Leave it to migrate for 15 minutes** before reading the result.

4/ READING AND INTERPRETATION OF THE RESULTS:

Wait **15 minutes for migration** before reading:

- **A NEGATIVE TEST:** a **single pink band** appears in the upper part of the strip (control band).

- **A POSITIVE TEST:** **two pink bands** appear on the strip (test band and control band).

The appearance of a test band after only 10 minutes of migration is considered a positive result. Any colouration of the test band, however slight, should be considered a positive result.

- **The absence of a control band in the upper part of the strip makes the test invalid.**

■ RECOMMENDATIONS

• STABILITY / STORAGE:

- 18 months at room temperature from the date of manufacture. The expiry date is indicated on the box and on the pouch of each testing device.
- Store at room temperature, between +2°C and +30°C. Avoid exposing the test to temperatures outwith this range.
- In case of refrigerated storage, let the reagents come to room temperature before use.

- It is recommended to keep the strips dry and to avoid touching the reactive parts of the membrane.

- **SAMPLES:**

- **The samples placed in the sampling tube cannot be stored for more than 2 hours.**
- When performing a test on a pooled sample, use a different sterile swab for each animal.

- **OTHER RECOMMENDATIONS:**

- During the test, the strip may show a light pink background coloration; this does not affect the quality of the results.
- **Never mix reagents from different batches.**

These recommendations are for guidance only; no diagnostic method is ever precise 100% of the time. The aim of this test is to detect BSRV antigens in order to help the veterinarian make an etiological diagnosis of respiratory disease in cattle. When interpreting the test, the veterinarian should always consider the medical history and clinical and epidemiological background of the animal. The final diagnosis is the prerogative and responsibility of the veterinarian performing the test.

Bio Véto Test and its distributors cannot be held responsible for the consequences of incorrect use or interpretation of the results given by this test.

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