**Speed Chlam™**

English	.....	1
Français	.....	6
Español	.....	11
Português	.....	16
Italiano	.....	21
Ελληνικά	.....	26
Deutsch	.....	31
Netherlands	.....	36
Procedure	.....	41

## FOR *IN VITRO* USE ONLY ENGLISH

### ■ CLINICAL APPLICATION

Chlamydiosis is a zoonosis that is particularly common in cats, birds, and farm animals. Chlamydiosis is caused by an intracellular Gram-negative bacterium, *Chlamydophila spp.*

Clinical signs are species-specific:

- Conjunctivitis, rhinitis, and pneumonia in cats.
- Diarrhea, respiratory insufficiency, and weight loss in birds.
- Abortion and pneumonia in cows, sheep, and goats.

Clinical signs are extremely variable and not specific; antigen testing is therefore essential to obtain a definitive diagnosis.

### ■ PRINCIPLE

Speed Chlam is a rapid qualitative immunoassay, using immunochromatography, for the detection of specific Chlamydial LPS (lipopolysaccharide) antigens. Samples (swabs) are taken from the conjunctiva, oropharynx, cloaca, or cervix, depending on the species.

Each sample is treated with extraction solution. The mixture is dropped on the test device. In the sample well, the stained particles of the conjugate bind to any Chlamydia antigens present in the sample. The resulting conjugate/antigen complexes migrate along the membrane via capillarity. They are captured by specific anti-Chlamydia antibodies that are bound to the membrane, forming a pink test band by accumulation of coloured particles. The mixture continues to migrate along the strip to the end of the membrane where the remaining coloured particles form a pink control band, which confirms the validity of the test.

## ■ PROCEDURE

### ► FOR EACH TEST YOU WILL NEED:

1 sterile swab, 1 extraction tube with filter, 1 test device, and the extraction solution vial.

The reagents should be used at room temperature.

**Do not mix reagents from different batches.**

### 1/ SAMPLE EXTRACTION:

- Take a swab from **conjunctival, oropharyngeal, or cloacal tissue**.
- Fill the extraction tube with **18 drops of extraction solution**.
- **Immerse the swab in the extraction tube**.
- **Swirl** the swab vigorously for **10 seconds**.
- Leave it to stand upright for **10 minutes at room temperature**.
- Remove the swab, pressing it against the wall of the extraction tube.

### 2/ ADD THE SAMPLE:

- Carefully fit the **filter dropper** to the vial.
- Hold the vial vertically and add **4 drops of extracted solution** to the sample well.

### 3/ READING AND INTERPRETATION OF RESULTS:

Wait **20 minutes** for the liquid to migrate before reading:

- **NEGATIVE TEST:** any test showing only **1 pink band** (control band) is a negative result.

- **POSITIVE TEST:** any test showing **2 distinct pink bands** (test band + control band) indicates a positive result.

The test is positive if a pink test band appears within 15 minutes of migration. Any colour change in the test band, even faint, should be considered as a positive result.

- **If no control band appears, the test is invalid.**

### ■ RECOMMENDATIONS

#### • SHELF LIFE:

- 24 months at room temperature (from manufacturing date). The expiry date is printed on each box and on each test pouch.
- Store at room temperature, between +2°C to +30°C. Keep the test away from sources of excessive heat or cold.

#### • SAMPLES:

- Use the provided swabs or synthetic swabs (rayon or dacron). Do not use swabs with cotton or calcium alginate tips, with wooden handles or impregnated with agar or gelatin.
- Remove excess mucous or blood from the endocervical area with a separate swab or cotton ball and discard. The sample must

contain as many Chlamydia infected cells as possible. Failure to remove mucus or blood may result in false positive results.

- If immediate testing is not possible, the sample can be placed in a dry plastic tube and stored in a refrigerator, between +2°C and +8°C, for up to 3 days; the test can detect both live and dead Chlamydia.

- **PRECAUTIONS:**

- To ensure correct migration of the sample, place the test kit on a horizontal surface.

- To ensure that the sample and reagent are correctly applied, always hold the extraction tube vertically.

- **Treat samples as a potential infectious biohazard.** Avoid contact with eyes or nose when collecting the samples. When the test is complete, discard used swabs in an appropriate hazardous waste container.

- **OTHER RECOMMENDATIONS:**

- During the test, the reading window may show a light pink background coloration; this does not affect the quality of the results.

- **Do not mix reagents from different batches.**

- If the kit is refrigerated, allow the reagents to come to room temperature before use.

- The extraction solution contains  $\text{NaN}_3$  which can become a toxic explosive acid on contact with lead or copper. Use large volumes of water to rinse reagents down the sink.

- **SAMPLING:**

- **Conjunctiva sampling:** lower the bottom eyelid and collect a sample with a swab. The sample must contain as many Chlamydia-infected conjunctival cells as possible.

- **Cervical sampling (bovine, ovine, caprine):** use special large animal swabs to remove excess mucus from the endocervical area and discard. Using one of the swabs supplied in the kit, rotate the swab against the tissues for 15 – 30 seconds, withdraw it in the transition zone to collect epithelial cells, without touching the vaginal epithelium.
- **Pharyngeal sampling:** insert the swab and rotate against the oropharyngeal epithelium to sample as many infected epithelial cells as possible.
- **Cloacal sampling:** use a swab to remove excess mucus from the cloaca and discard. Insert another swab (provided in the kit) into the cloaca and rotate to sample infected epithelial cells.
- **Bird dropping:** plunge a swab into the droppings (dropping can be stored in a bag or a box).
- **Organ sampling (birds) and placental tissues (bovine, ovine, and caprine):** simply swab the surface of the tissue (liver, lung, placenta, etc.).

The above recommendations are only guidelines; no test is 100% accurate at all times and under all conditions. The purpose of this test is to assist practicing vets with the diagnosis of Chlamydiosis by the detection of *Chlamydia psittaci* antigens. All test results should be interpreted in the light of the patient's history, clinical examination, and the results of any further diagnostic tests. The definitive diagnosis remains the veterinarian's prerogative and responsibility.

Bio Veto Test and its distributors cannot be held responsible for the consequences of misuse or misinterpretation of the results of the test.

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