**Speed Parvo™**

# Speed Parvo™

## Veterinary diagnostic kit

For *in vitro* use only

### ■ CLINICAL APPLICATION

Canine Parvovirus causes acute, highly infectious, and severely debilitating disease, especially in puppies under 6 months of age. Symptoms include haemorrhagic gastroenteritis associated with panleucopenia, anorexia, dehydration, and depression. However, these symptoms are not specific to canine parvovirus infection thus clinical signs alone cannot be used to differentiate parvovirosis from other intestinal diseases in dogs. Speed Parvo provides the veterinarian with a diagnosis of canine parvovirus only 2 or 3 days after infection.

Maternal antibodies protect puppies during their first 5 to 6 weeks of life, during which time vaccination is inactive. The puppies' immune system only starts to mature from 8 weeks of age, so between 6 and 13 weeks puppies are highly susceptible to infection.

The virus is shed in the faeces of infected animals and remains actively infectious in the environment for up to 8 months. Puppies can be infected via the oronasal route, which is why some kennels are permanent carriers of the disease. Early diagnosis of canine parvovirus infection enables the veterinarian to administer appropriate life-saving treatment and to quarantine infected dogs.

### ■ PRINCIPLE

The Speed Parvo test is a rapid qualitative immunoassay, using immunochromatography, for the detection of canine parvovirus antigens.

For each test, in the sample well, the stained particles of the conjugate bind to any parvovirus antigens present in the sample. The resulting conjugate/antigen complexes migrate along the membrane via capillarity. They are captured by specific antibodies that are bound to the membrane, forming a pink test band by accumulation of coloured particles. The mixture continues to migrate along the strip to the end of the membrane where the remaining coloured particles form a pink control band, which confirms the validity of the test.

## ■ PROCEDURE

### ► FOR EACH TEST YOU WILL NEED:

1 test device, 1 swab, 1 filter cap, and 1 reagent vial

The reagents should be used at room temperature.

**Do not mix reagents from different batches.**

### 1/ ADD THE SAMPLE:

- Coat the swab with faeces from a rectal or faecal sample. Discharge the swab in the reagent vial.
- Carefully fit the filter dropper to the vial.
- Shake the vial thoroughly to break up and homogenise the particles.

### 2/ PERFORMING THE TEST:

Hold the vial vertically and add **6 drops** to the sample well.

### 3/ READING AND INTERPRETATION:

**Wait 10 minutes** for the liquid to migrate before reading:

- **NEGATIVE TEST:** any test showing only **1 pink band** (control band) is a negative result.

- **POSITIVE TEST:** any test showing **2 distinct pink bands** (test band + control band) indicates a positive result. Any coloration of the test band (even a faint red coloration) must be considered as a positive result.

- **If no control band appears, the test is invalid.**

### ■ RECOMMENDATIONS

#### • SHELF LIFE:

- 24 months at room temperature (from manufacturing date). The expiry date is printed on each box and on each test pouch.
- Store at room temperature between +2°C to +30°C. Keep the test away from excess heat or cold.

#### • SAMPLES:

- Speed Parvo is only suitable for use with faecal or rectal swabs.
- Faecal samples can be stored for up to 14 days at between +2°C and +8°C.
- If they need to be stored for longer, the contaminated reagents vial can be stored for several months at -20°C.

#### • PRECAUTIONS:

- To ensure correct migration of the sample, place the test device on a horizontal surface.
- To ensure that the sample is applied correctly, always hold the

reagent bottle vertically.

• **OTHER RECOMMENDATIONS:**

- The reading window may show a light pink background coloration; this will not affect the accuracy of the results.
- **Do not mix reagents from different batches.**
- Allow the reagents to come to room temperature before use.
- Performance of a test within 3 weeks after vaccination can sometimes lead to the appearance of a weak positive reaction due to the transient presence of vaccination antigens in the faeces.
- All instruments and materials used for sampling should be considered as potentially infectious and be disposed of in the hazardous waste container.

The above recommendations are only guidelines; no test is 100% accurate at all times and under all conditions. The purpose of this test is to detect Parvovirus antigens to diagnose parvovirus. All test results should be interpreted in the light of the patient's history, clinical examination, and the results of any further diagnostic tests. The definitive diagnosis remains the veterinarian's prerogative and responsibility.

Bio Veto Test and its distributors cannot be held responsible for the consequences of misuse or misinterpretation of the results given by the test.

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