# *Speed Duo*™ Leish K/Ehrli

## Speed Duo™ Leish K/Ehrli

Veterinary diagnostic kit  
For *in vitro* use only

### ■ CLINICAL APPLICATION

Canine leishmaniosis is caused by the proliferation of a parasite called *Leishmania infantum*, in the macrophages of its host. The disease is transmitted by a little spinning dipter: the phlebotome or sand fly.

Canine ehrlichiosis is a bacterial infection caused by a *Rickettsiae*, *Ehrlichia canis*, transmitted by a tick: *Rhipicephalus sanguineus*.

These vector-borne diseases are found worldwide. Their geographical distribution is related to that of their vectors. In Europe, leishmaniosis and ehrlichiosis are very common around the Mediterranean basin and co-infections are frequently diagnosed.

Leishmaniosis and ehrlichiosis present many clinical similarities and the differential diagnosis can be challenging:

- Following infection, symptoms often go unnoticed by the owner.
- Asymptomatic or sub-clinical phases are common and often last a long time.
- Sick dogs present very similar clinical symptoms: depression, weight loss, weakness, lymph node enlargement...

The possibility of concurrent infection and the clinical similarity between the two diseases makes the serological diagnosis of leishmaniosis and ehrlichiosis very important. In endemic areas, regular screening of the two diseases is strongly recommended to enable early detection in asymptomatic dogs.

EN-1

## ■ PRINCIPLE

Speed Duo Leish K/Ehrli is a rapid qualitative immunoassay, using immunochromatography, for the detection of anti-*Leishmania infantum* and anti-*Ehrlichia canis* antibodies in dogs.

For each test, one drop of the sample (serum, plasma or anti-coagulated whole blood) is added to the sample wells. In the sample well, the stained particles of the conjugate bind to any anti-*Leishmania infantum* antibodies and anti-*Ehrlichia canis* antibodies present in the sample. The resulting conjugate/anti-*Leishmania infantum* antibody and conjugate/anti-*Ehrlichia canis* antibody complexes migrate along the membrane via capillarity. They are captured by specific *Leishmania* antigens (Leish part) and specific *Ehrlichia* antigens (Ehrli part) that are bound to the membrane, forming a pink test band by accumulation of coloured particles. The mixture continues to migrate along the strip to the end of the membrane where the remaining coloured particles form a pink control band, which confirm the validity of the test.

## ■ PROCEDURE

### ► FOR EACH TEST YOU WILL NEED:

1 test device, 1 single-use pipette and the reagent dropper bottle.  
The reagents should be used at room temperature

**Do not mix reagents from different batches**

### 1/ ADD THE SAMPLE:

Serum, plasma, anti-coagulated whole blood (EDTA or heparin):  
Using the provided single-use pipette held vertically, **add 1 drop** of sample **at the centre of each sample well and wait 10 to 15 seconds.**

### 2/ ADD THE REAGENT:

- **Add 5 drops of reagent** to each sample well, holding the dropper bottle in a vertical position.
- If no migration is observed within 2 minutes, **add an extra 2 drops** of reagent to the sample well.

EN-2

### 3/ READING AND INTERPRETATION OF RESULTS:

**Wait 15 minutes** for migration and read:

- **NEGATIVE TEST:** any test showing only **1 pink band** (control band) is a negative result.

- **POSITIVE TEST:** any test showing **2 distinct pink bands** (test band + control band) indicates a positive result. The test is positive if a pink test band appears within 10 minutes of migration. Any colour change in the test band, even faint, should be considered as a positive result.

- **If no control band appears, the test is invalid.**

### ■ RECOMMENDATIONS

#### • SHELF LIFE:

- 24 months at room temperature (from manufacturing date). The expiry date is printed on each box and on each test pouch.
- Store at room temperature, between +2°C to +30°C. Keep the test away from sources of excessive heat or cold.

#### • SAMPLES:

- Serum, plasma or anti-coagulated whole blood has to be used at room temperature.
- For testing within 24 hours of sampling: store serum, plasma or anti-coagulated whole blood at room temperature.
- For testing between 24 to 72 hours of sampling: store serum or plasma between +2°C and +8°C (refrigerator). Whole blood with anticoagulant storage between +2°C and +8°C for several hours may induce haemolysis, (not visible to the naked eye), which may interfere with the interpretation of the results.

EN-3

- Testing more than 72 hours after sampling: store serum or plasma at -20°C (freezer). Never freeze a whole blood sample.

• **PRECAUTIONS:**

- To ensure correct migration of the sample, place the test device on a horizontal surface.
- To standardize the amount of sample recommended for each test, use only the single-use pipettes supplied in the kit.
- To ensure that the sample and reagent are correctly applied, always hold the pipette or dropper bottle vertically.

• **OTHER RECOMMENDATIONS:**

- The reading window may show a light pink background coloration; this does not affect the quality of the results.
- **Do not mix reagents from different batches.**
- Allow the reagents to come to room temperature before use.
- Use a new single-use pipette for each test.

The above recommendations are only guidelines; no test is 100% accurate at all times and under all conditions. The purpose of this test is to detect anti-*Leishmania infantum* and anti-*Ehrlichia canis* antibodies to diagnose leishmaniosis and ehrlichiosis respectively. All test results should be interpreted in the light of patient's history, clinical examination and the results of any further diagnostic tests. The definitive diagnosis remains the veterinarian's prerogative and responsibility.

Bio Veto Test and its distributors cannot be held responsible for the consequences of misuse or misinterpretation of the results given by the test.

**Kit diagnostic vétérinaire**  
**Usage *in vitro* uniquement****■ INTERET CLINIQUE**

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**SPEED DUO, *Speed* DUO, SPEED and *Speed* are registered trademarks or trademarks of Virbac or its affiliates**
SPEED DUO, *Speed* DUO, SPEED & *Speed* zijn geregistreerde trademarks of trademarks van Virbac en zijn filialen