*Speed*<sup>TM</sup> FeLV

**Veterinary diagnostic kit  
For *in vitro* use only****■ CLINICAL APPLICATION**

Feline Leukaemia virus (FeLV) is a retrovirus that is present worldwide. It is primarily transmitted through close contact with infected cats via saliva, nasal secretions, or faeces.

Following infection, the initial viraemia is handled in one of two ways:

- the immune system controls the infection and quickly eliminates the circulating virus.
- the viraemia persists for several months or even years without any clinical symptoms. Such animals are considered to be chronic carriers and will excrete the virus intermittently, thus providing a source of infection for other cats.

The clinical signs of FeLV infection are polymorphic and non-specific and include depression, anorexia, pyrexia, immunodeficiency syndrome, and anaemia, etc. Further diagnostic tests are required to obtain a definitive diagnosis.

Regular FeLV testing can be used to determine disease status and to indicate whether specific precautions need to be taken (vaccination of healthy cats, neutering, and medical supervision of infected or sick cats...).

**■ PRINCIPLE**

Speed FeLV is a rapid qualitative immunoassay, using immunochromatography, for the detection of FeLV P27 antigen in cats.

For each test, one drop of the sample (serum, plasma, or anticoagulated whole blood) is added to the sample well. In the sample well, the stained particles of the conjugate bind to any FeLV P27 antigens present in the sample. The resulting conjugate/antigen complexes migrate along the membrane via capillarity. They are captured by specific anti-FeLV antibodies that are bound to the membrane, forming a pink test band by accumulation of coloured particles. The mixture continues to migrate along the strip to the end of the membrane where the remaining coloured particles form a pink control band, which confirms the validity of the test.

## ■ PROCEDURE

### ► FOR EACH TEST YOU WILL NEED:

1 test device, 1 single-use pipette, and the reagent dropper bottle.

The reagents should be used at room temperature.

**Do not mix reagents from different batches.**

#### 1/ ADD THE SAMPLE:

Serum, plasma, anti-coagulated whole blood (EDTA or heparin); using the provided single-use pipette held vertically, **add 1 drop of sample at the centre of the sample well and wait 10 to 15 seconds.**

#### 2/ ADD THE REAGENT:

- **Add 5 drops of reagent** to the sample well, holding the dropper bottle in a vertical position.
- If no migration is observed within 2 minutes, **add an extra 2 drops** of reagent to the sample well.

### **3/ READING AND INTERPRETATION OF RESULTS:**

**Wait 15 minutes** for migration and read:

- **NEGATIVE TEST:** any test showing only **1 pink band** (control band) is a negative result.

- **POSITIVE TEST:** any test showing **2 distinct pink bands** (test band + control band) indicates a positive result. The test is positive if a pink test band appears within 10 minutes of migration. Any colour change in the test band, even faint, should be considered as a positive result.

- **If no control band appears, the test is invalid.**

### **■ RECOMMENDATIONS**

#### **• SHELF LIFE:**

- 24 months at room temperature (from manufacturing date). The expiry date is printed on each box and on each test pouch.
- Store at room temperature, between + 2°C to + 30°C. Keep the test away from sources of excessive heat or cold.

#### **• SAMPLES:**

- Serum, plasma, or anticoagulated whole blood must be used at room temperature.
- Testing within 24 hours of sampling: store serum, plasma, or anticoagulated whole blood at room temperature.
- Testing between 24 to 72 hours of sampling: store serum or plasma between + 2°C to + 8°C (refrigerator). Storing whole blood between + 2°C to + 8°C for several hours may induce haemolysis (not visible to the naked eye), which may interfere with the interpretation of the results.

- Testing more than 72 hours after sampling; store serum or plasma at -20°C (freezer). Never freeze a whole blood sample.

- **PRECAUTIONS:**

- To ensure correct migration of the sample, place the test device on a horizontal surface.
- To standardize the amount of sample recommended for each test, only use the single-use pipettes supplied in the kit.
- To ensure that the sample and reagent are correctly applied, always hold the pipette or the dropper bottle vertically.

- **OTHER RECOMMENDATIONS:**

- The reading window may show a light pink background colouration; this does not affect the quality of the results.
- **Do not mix reagents from different batches.**
- Allow the reagents to come to room temperature before use.
- Use a new single-use pipette for each test.

The above recommendations are only guidelines; no test is 100% accurate at all times and under all conditions. The purpose of this test is to detect FeLV antigens to aid in the diagnosis of FeLV. All test results should be interpreted in the light of the patient's history, clinical examination, and the results of any further diagnostic tests. The definitive diagnosis remains the veterinarian's prerogative and responsibility.

Bio Veto Test and its distributors cannot be held responsible for the consequences of misuse or misinterpretation of the results of the test.

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SPEED and *Speed* are trademarks of Virbac or its affiliates