Speed[™] Parvo

Speed™ Parvo

Veterinary diagnostic kit For *in vitro* use only

■ CLINICAL APPLICATION

Canine Parvovirus causes acute, highly infectious, and severely debilitating disease, especially in puppies under 6 months of age. Symptoms include haemorrhagic gastroenteritis associated with panleucopenia, anorexia, dehydration, and depression. However, these symptoms are not specific to canine parvovirus infection thus clinical signs alone cannot be used to differentiate parvovirus from other intestinal diseases in dogs. Speed Parvo provides the veterinarian with a diagnosis of canine parvovirus only 2 or 3 days after infection.

Maternal antibodies protect puppies during their first 5 to 6 weeks of life, during which time vaccination is not recommended due to maternal antibody interference. The puppies' immune system only starts to mature from 8 weeks of age, so between 6 and 16 weeks puppies are highly susceptible to infection.

The virus is shed in the faeces of infected animals and remains actively infectious in the environment for up to 8 months. Puppies can be infected via the oronasal route, which is why some kennels are permanent carriers of the disease. Early diagnosis of canine parvovirus infection enables the veterinarian to administer appropriate life-saving treatment and to quarantine infected dogs.

■ PRINCIPLE

The Speed Parvo test is a rapid qualitative immunoassay, using immunochromatography, for the detection of canine parvovirus antigens.

For each test, in the sample well, the stained particles of the conjugate

bind to any parvovirus antigens present in the sample. The resulting conjugate/antigen complexes migrate along the membrane via capillarity. They are captured by specific antibodies that are bound to the membrane, forming a pink test band by accumulation of coloured particles. The mixture continues to migrate along the strip to the end of the membrane where the remaining coloured particles form a pink control band, which confirms the validity of the test.

■ PROCEDURE

► FOR EACH TEST YOU WILL NEED:

1 test device, 1 swab and 1 reagent vial.

The reagents should be used at room temperature.

Do not mix reagents from different batches.

1/ ADD THE SAMPLE:

Faecal samples should be collected only with the provided swabs. Do not use the top of the cap of the reagent vial to collect the samples.

- Coat **entirely** the fibrous tip of the swab with faeces from a rectal or faecal sample and discharge it into the reagent vial.
- Shake the swab vigorously in the vial for a few seconds and then remove the maximum amount of faeces by pressing and turning the fibrous tip of the swab against the walls of the vial.
- Close the vial and shake well to homogenize the content.

2/ PERFORMING THE TEST:

- Break the detachable top of the cap to release the inverted dropper tip.
- Invert the vial and add **5 drops** of the content to the sample well.
- Leave to migrate.

EN-2

3/ READING AND INTERPRETATION:

Wait **10 minutes** for the liquid to migrate before reading:

- **NEGATIVE TEST:** any test showing only **1 pink band** (control band) is a negative result.

- **POSITIVE TEST:** any test showing **2 distinct pink bands** (test band + control band) indicates a positive result. Any coloration of the test band (even a faint red coloration) must be considered as a positive result.

- **If no control band appears, the test is invalid.**

■ RECOMMENDATIONS

• SHELF LIFE:

- 24 months at room temperature (from manufacturing date). The expiry date is printed on each box and on each test pouch.
- Store at room temperature between +2°C to +30°C. Keep the test away from excess heat or cold.

• SAMPLES:

- Speed Parvo is only suitable for use with faecal or rectal swabs.
- It is recommended to **entirely** coat the fibrous tip of the swab with a thin layer of faeces
- Faecal samples can be stored for up to 14 days at between +2°C and +8°C.
- If they need to be stored for longer, the contaminated reagent's vials can be stored for several months at -20°C.

EN-3

- **PRECAUTIONS:**

- To ensure correct migration of the sample, place the test device on a horizontal surface.
- To ensure that the sample is applied correctly, always hold the reagent bottle vertically.
- It is advisable to wear protective gloves, protective clothing when taking the sample and during all stages of the test procedure.
- Perform the test procedures on a dry and clean work surface.
- The fecal specimen and all the materials used for the test procedure should be considered as potentially infectious and they should be disposed of in accordance with local regulations.

- **OTHER RECOMMENDATIONS:**

- The reading window may show a light pink background coloration; this will not affect the accuracy of the results.
- **Do not mix reagents from different batches.**
- Allow the reagents to come to room temperature before use.
- Performance of a test within 3 weeks after vaccination can sometimes lead to the appearance of a weak positive reaction due to the transient presence of vaccination antigens in the faeces.

The above recommendations are only guidelines; no test is 100% accurate at all times and under all conditions. The purpose of this test is to detect Parvovirus antigens to diagnose parvovirus. All test results should be interpreted in the light of the patient's history, clinical examination, and the results of any further diagnostic tests. The definitive diagnosis remains the veterinarian's prerogative and responsibility.

Bio Veto Test and its distributors cannot be held responsible for the consequences of misuse or misinterpretation of the results given by the test.

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SPEED and *Speed* are trademarks of Virbac or its affiliates