
Speed™ Progesterone

Veterinary diagnostic test For *in vitro* use only


Fluorescence immunoassay test for the quantitative measurement of progesterone in canine serum or plasma.

1


2

100µL


3


4

X5


5


6

2 X

7

1

2

CLINICAL APPLICATION
The exact timing of ovulation in bitches can be highly variable: 20% of bitches ovulate outside of the 10th to 14th day of heat. Vaginal smears, combined with progesterone assays, provide the most accurate determination of ovulation to enable scheduling of insemination or mating. Progesterone level is low prior to ovulation. It begins to increase at LH peak, reaching 5 to 10 ng/mL at the time of ovulation.

PRINCIPLE
The Speed™ Progesterone assay is a rapid immunochromatographic quantitative test which measures the circulating progesterone by laser-induced fluorescence analysis.

PRECAUTIONS AND WARNINGS
For use by veterinarians only

SAMPLE

- Use only serum or plasma (heparin-Li) samples. Do not use whole blood or EDTA treated samples.
- Do not use tubes containing serum separator gel.
- No interference with haemolysis has been seen.
- Highly lipaemic samples may affect the results. A 6 hour fasting prior to sampling is recommended to limit hyperlipaemia.
- If the test is not performed within an hour, transfer the serum or plasma to a plain tube.
- Samples can be stored in a refrigerator (+2°C to +8°C) for up to 72 hours.
- For testing beyond 72 hours, the serum or plasma should be stored in a freezer (-20°C).
- Stored samples should be allowed to reach room temperature (18°C to 27°C) and be centrifuged before analysis.
- The presence of fibrin strands in the sample may lead to erroneous assay results due to pipette clogging and reduced volume of aspirated sample. In this case, recentrifuge the sample.

PROPER PIPETTE USAGE

- Hold the pipette, the sample and reagent tubes vertically at all times, to ensure the correct sampling volume.
- Use a new disposable pipette tip for each test.
- Press and hold the pipette plunger and insert the tip into the sample, being careful not to touch the walls of the tube. Release the pipette plunger slowly while keeping the end of the tip in the sample.
- Remove the pipette from the sample tube and check for air bubbles in the tip. If air bubbles are present, discard the sample and repeat the sampling.

REAGENT AND SAMPLE MIXTURE

- To ensure correct migration of the sample, place the test device on a horizontal surface.
- The reagent tube contains a standard quantity of reagent. Any alteration of this quantity can lead to erroneous assay results.
- Use the sample/reagent mixture immediately after reconstitution.
- Do not keep the mixture for subsequent analysis.

RECOMMENDATIONS

- Before using the first test of each box, register the lot calibration chip in the Speed Reader™ analyser
- The expiry date is indicated on the box and on each test pouch.
- Store the tests and the reagents between +2°C and +8°C.
- The specimen and all the materials used for the test procedure should be considered as potentially infectious and they should be disposed of in accordance with local regulations.
- Use appropriate and clean protective equipment (gloves and gown)
- In case of skin or eye contact with the reagent, rinse immediately with clean water and ask for medical advice.

TEST PROCEDURE

The Speed™ Progesterone test cannot be visually interpreted and should be read only with a Speed Reader™ analyser. For complete instructions, please refer to the Speed Reader™ analyser operator's guide.

usage. Be careful not to touch the walls of the reagent tube with the tip

- **Carefully seal the reagent tube with the dropper cap**
- **Gently mix the contents by inverting the tube at least 5 times.** DO NOT SHAKE

SAMPLE APPLICATION

- Remove the plastic cover from the dropper cap and **discard the first 2 drops** of the mixture, to wash out the dead space of the dropper cap
- Completely remove the cartridge from Speed Reader™ analyser and add **2 drops** of the mixture into the sample well of the test device, by holding the vial in an **absolutely vertical position**.

READING

- Insert the test device into the cartridge holder, when prompted by the analyser, in the direction indicated by the arrow (the sample well first) and follow the on-screen instructions.
- Read the result on the screen

INTERPRETATION
Progesterone concentration is measured in ng/mL and nmol/L. - Dynamic Range: 1 ng/mL to 20 ng/mL and 3.18 nmol/L to 63.60 nmol/L

Concentration	Interpretation	Action**
0 to 2 ng/mL 0 to 6,36 nmol/L	Baseline progesterone	Repeat the test every 2 days until the progesterone is ≥ 2 ng/mL ($\geq 6,36$ nmol/L)
2 to 5 ng/mL 6,36 to 15,90 nmol/L	LH peak. Progesterone level is increasing	Repeat the test every day until the progesterone is ≥ 5 ng/mL ($\geq 15,90$ nmol/L)
5 to 10 ng/mL 15,90 to 31,80 nmol/L	Consistent with ovulation period	Mating or insemination is recommended twice between 1 and 4 days later.
> 10 ng/mL >31,80 nmol/L	The fertile period might be close to the end	Mate or inseminate immediately. Caution: the fertile period might be over

Conversion: 1 ng/mL = 3,18 nmol/L
* Reference values represent mean progesterone levels during the oestrus cycle, considerable variations within normal may often occur.
** The veterinarian should interpret all test results in light of the patient's history, ultrasonography, and/or vaginal smears.

These recommendations are for guidance only. The aim of this test is to measure canine progesterone blood concentrations. Interpretation of the result by the veterinarian should always take into account the history, clinical examination, and any further diagnostic test results, as no diagnostic method is 100% accurate. The definitive diagnosis is the prerogative and responsibility of the veterinarian.
Bio Veto Test and its distributors cannot be held responsible for any consequences linked to incorrect use of this test or misinterpretation of the results.

SPEED READER, *Speed Reader*, SPEED and *Speed* are registered trademarks or trademarks of Virbac or its affiliates

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